

OR211P Antibody



PACO04358

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR211P

Uniprot

Q8NGU4

Synonyms:

OR211P; OR212; OR213P; OR214P; Putative olfactory receptor 211; Putative olfactory receptor 212; Putative olfactory receptor 213; Putative olfactory receptor 214

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 211.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A