

GPR143 Antibody



PACO04476

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

GPR143

Uniprot

P51810

Synonyms:

GPR143; OA1; G-protein coupled receptor 143; Ocular albinism type 1 protein

Immunogen:

Synthesized peptide derived from the Internal region of human GPR143.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A