

TAS2R3 Antibody



PACO04505

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

TAS2R3

Uniprot

Q9NYW6

Synonyms:

TAS2R3; Taste receptor type 2 member 3; T2R3

Immunogen:

Synthesized peptide derived from the Internal region of human T2R3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A