

# TAS2R13 Antibody



PACO04509

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

TAS2R13

**Uniprot**

Q9NYV9

**Synonyms:**

TAS2R13; Taste receptor type 2 member 13; T2R13; Taste receptor family B member 3; TRB3

**Immunogen:**

Synthesized peptide derived from the Internal region of human T2R13.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A