

TAS2R10 Antibody



PACO04510

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

TAS2R10

Uniprot

Q9NYW0

Synonyms:

TAS2R10; Taste receptor type 2 member 10; T2R10; Taste receptor family B member 2; TRB2

Immunogen:

Synthesized peptide derived from the Internal region of human T2R10.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A