

TAS2R14 Antibody



PACO04511

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

TAS2R14

Uniprot

Q9NYV8

Synonyms:

TAS2R14; Taste receptor type 2 member 14; T2R14; Taste receptor family B member 1; TRB1

Immunogen:

Synthesized peptide derived from the Internal region of human T2R14.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A