

# CYP39A1 Antibody



PACO04532

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:****Gene ID:**

CYP39A1

**Uniprot**

Q9NYL5

**Synonyms:**

CYP39A1; 24-hydroxycholesterol 7-alpha-hydroxylase; Cytochrome P450 39A1;  
hCYP39A1; Oxysterol 7-alpha-hydroxylase

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human CYP39A1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A