

SPTBN5 Antibody



PACO04533

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300, IHC:1:200-1:1000

Protein Background:

Gene ID:

SPTBN5

Uniprot

Q9NRC6

Synonyms:

SPTBN5; BSPECV; HUBSPECV; HUSPECV; Spectrin beta chain; non-erythrocytic 5; Beta-V spectrin

Immunogen:

Synthesized peptide derived from the Internal region of human Spectrin Beta V.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A