

VPS54 Antibody



PACO04541

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

VPS54

Uniprot

Q9P1Q0

Synonyms:

VPS54; HCC8; Vacuolar protein sorting-associated protein 54; Hepatocellular carcinoma protein 8; Tumor antigen HOM-HCC-8; Tumor antigen SLP-8p

Immunogen:

Synthesized peptide derived from the Internal region of human VPS54.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A