

PIGF Antibody



PACO04586

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

PIGF

Uniprot

Q07326

Synonyms:

PIGF; Phosphatidylinositol-glycan biosynthesis class F protein; PIG-F; GPI11 homolog

Immunogen:

Synthesized peptide derived from the C-terminal region of human PIG-F.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A