

# POLE Antibody



PACO04632

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

POLE

**Uniprot**

Q07864

**Synonyms:**

POLE; POLE1; DNA polymerase epsilon catalytic subunit A; DNA polymerase II subunit A

**Immunogen:**

Synthesized peptide derived from the Internal region of human DNA pol epsilon A.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A