

GPR85 Antibody



PACO04634

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

GPR85

Uniprot

P60893

Synonyms:

GPR85; SREB2; Probable G-protein coupled receptor 85; Super conserved receptor expressed in brain 2

Immunogen:

Synthesized peptide derived from the Internal region of human GPR85.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A