

TAS2R5 Antibody



PACO04639

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

TAS2R5

Uniprot

Q9NYW4

Synonyms:

TAS2R5; Taste receptor type 2 member 5; T2R5

Immunogen:

Synthesized peptide derived from the Internal region of human T2R5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A