

SH3TC1 Antibody



PACO04640

Product Information

Size: 50ug	Protein Background:
Reactivity: Human	Gene ID: SH3TC1
Source: Rabbit	Uniprot Q8TE82
Isotype: IgG	Synonyms: SH3TC1; SH3 domain and tetratricopeptide repeat-containing protein 1
Applications: ELISA, WB	Immunogen: Synthesized peptide derived from the Internal region of human SH3TC1.
Recommended dilutions: WB:1:500-1:2000	Storage: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A