

FAM35A Antibody



PACO04646

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

FAM35A

Uniprot

Q86V20/A6NDK8

Synonyms:

FAM35A; Protein FAM35A; FAM35B; Protein FAM35B

Immunogen:

Synthesized peptide derived from the Internal region of human FAM35A/B.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A