

RNF111 Antibody



PACO04658

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

RNF111

Uniprot

Q6ZNA4

Synonyms:

RNF111; E3 ubiquitin-protein ligase Arkadia; RING finger protein 111

Immunogen:

Synthesized peptide derived from the C-terminal region of human Arkadia.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A