

ZNF280D Antibody



PACO04659

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

ZNF280D

Uniprot

Q6N043

Synonyms:

ZNF280D; KIAA1584; SUHW4; ZNF634; Zinc finger protein 280D; Suppressor of hairy wing homolog 4; Zinc finger protein 634

Immunogen:

Synthesized peptide derived from the Internal region of human SUHW4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A