

# ARHGAP17 Antibody



PACO04685

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

ARHGAP17

**Uniprot**

Q68EM7

**Synonyms:**

ARHGAP17; RICH1; MSTP066; MSTP110; Rho GTPase-activating protein 17; Rho-type GTPase-activating protein 17; RhoGAP interacting with CIP4 homologs protein 1; RICH-1

**Immunogen:**

Synthesized peptide derived from the Internal region of human ARHGAP17.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A