

DNAJC17 Antibody



PACO04689

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

DNAJC17

Uniprot

Q9NVM6

Synonyms:

DNAJC17; DnaJ homolog subfamily C member 17

Immunogen:

Synthesized peptide derived from the N-terminal region of human DnaJ C17.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A