

# CDK5RAP2 Antibody



PACO04720

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

CDK5RAP2

**Uniprot**

Q96SN8

**Synonyms:**

CDK5RAP2; CEP215; KIAA1633; CDK5 regulatory subunit-associated protein 2; CDK5 activator-binding protein C48; Centrosome-associated protein 215

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human CDK5 Activator-binding C48.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A