

# PRKCD Antibody



PACO04732

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

PRKCD

**Uniprot**

Q05655

**Synonyms:**

PRKCD; Protein kinase C delta type; Tyrosine-protein kinase PRKCD; nPKC-delta

**Immunogen:**

Synthesized peptide derived from human PKC delta around the non-phosphorylation site of Y52.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A