

# Phospho-PRKCE (S729) Antibody



PACO04735

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

PRKCE

**Uniprot**

Q02156

**Synonyms:**

PRKCE; PKCE; Protein kinase C epsilon type; nPKC-epsilon

**Immunogen:**

Synthesized peptide derived from human PKC epsilon around the phosphorylation site of S729.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A