

PACO04789

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

NCLN

Uniprot

Q969V3

Synonyms:

NCLN; Nicalin; Nicastrin-like protein

Immunogen:

Synthesized peptide derived from the C-terminal region of human Nicalin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A