

SHD Antibody



PACO04791

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

SHD

Uniprot

Q96IW2

Synonyms:

SHD; SH2 domain-containing adapter protein D

Immunogen:

Synthesized peptide derived from the N-terminal region of human Shd.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A