

# DDX24 Antibody



PACO04793

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

DDX24

**Uniprot**

Q9GZR7

**Synonyms:**

DDX24; ATP-dependent RNA helicase DDX24; DEAD box protein 24

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human DDX24.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A