

INTS2 Antibody



PACO04821

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

INTS2

Uniprot

Q9H0H0

Synonyms:

INTS2; KIAA1287; Integrator complex subunit 2; Int2

Immunogen:

Synthesized peptide derived from the Internal region of human INTS2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A