

Phospho-ARHGAP31 (T789) Antibody



PACO04822

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

ARHGAP31

Uniprot

Q9ULL6

Synonyms:

ARHGAP31; CDGAP; KIAA1204; Rho GTPase-activating protein 31; Cdc42 GTPase-activating protein

Immunogen:

Synthesized peptide derived from human CdGAP around the phosphorylation site of T789.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A