

COL20A1 Antibody



PACO04833

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

COL20A1

Uniprot

Q9P218

Synonyms:

COL20A1; KIAA1510; Collagen alpha-1(XX) chain

Immunogen:

Synthesized peptide derived from the C-terminal region of human COL20A1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A