

RAD23A Antibody



PACO04865

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

RAD23A

Uniprot

P54725

Synonyms:

RAD23A; UV excision repair protein RAD23 homolog A; HR23A; hHR23A

Immunogen:

Synthesized peptide derived from the Internal region of human Rad23A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A