

Phospho-RB1 (S249) Antibody



PACO04889

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:**Gene ID:**

RB1

Uniprot

P06400

Synonyms:

RB1; Retinoblastoma-associated protein; p105-Rb; pRb; Rb; pp110

Immunogen:

Synthesized peptide derived from human Rb around the phosphorylation site of S249.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A