

RELB Antibody



PACO04907

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

RELB

Uniprot

Q01201

Synonyms:

RELB; Transcription factor RelB; I-Rel

Immunogen:

Synthesized peptide derived from human RelB around the non-phosphorylation site of S552.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A