

# BCL7A Antibody



PACO04927

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

BCL7A

**Uniprot**

Q4VC05

**Synonyms:**

BCL7A; B-cell CLL/lymphoma 7 protein family member A

**Immunogen:**

Synthesized peptide derived from the Internal region of human Bcl-7a.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A