

BPIFA4P Antibody



PACO04931

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

BPIFA4P

Uniprot

Q86YQ2

Synonyms:

LATH; BASE; Putative latherin; Breast cancer and salivary gland-expressed protein

Immunogen:

Synthesized peptide derived from the Internal region of human Latherin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A