

RPL3L Antibody



PACO04937

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

RPL3L

Uniprot

Q92901

Synonyms:

RPL3L; 60S ribosomal protein L3-like

Immunogen:

Synthesized peptide derived from the Internal region of human Ribosomal Protein L3L.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A