

RPL37 Antibody



PACO04945

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

RPL37

Uniprot

P61927

Synonyms:

RPL37; 60S ribosomal protein L37; G1.16

Immunogen:

Synthesized peptide derived from the Internal region of human Ribosomal Protein L37.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A