

RPS15 Antibody



PACO04967

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

RPS15

Uniprot

P62841

Synonyms:

RPS15; RIG; 40S ribosomal protein S15; RIG protein

Immunogen:

Synthesized peptide derived from the Internal region of human Ribosomal Protein S15.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A