

RPS21 Antibody



PACO04971

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

RPS21

Uniprot

P63220

Synonyms:

RPS21; 40S ribosomal protein S21

Immunogen:

Synthesized peptide derived from the Internal region of human Ribosomal Protein S21.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A