

RPS25 Antibody



PACO04974

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

RPS25

Uniprot

P62851

Synonyms:

RPS25; 40S ribosomal protein S25

Immunogen:

Synthesized peptide derived from the C-terminal region of human Ribosomal Protein S25.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A