

# S100A5 Antibody



PACO04980

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

S100A5

**Uniprot**

P33763

**Synonyms:**

S100A5; S100D; Protein S100-A5; Protein S-100D; S100 calcium-binding protein A5

**Immunogen:**

Synthesized peptide derived from the Internal region of human S-100A5.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A