

SH2D4A Antibody



PACO05005

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:

Gene ID:

SH2D4A

Uniprot

Q9H788

Synonyms:

SH2D4A; SH2A; SH2 domain-containing protein 4A; Protein SH(2)A

Immunogen:

Synthesized peptide derived from the Internal region of human SH2D4A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A