

CIDEC Antibody



PACO05006

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

CIDEC

Uniprot

Q96AQ7

Synonyms:

CIDEC; FSP27; Cell death activator CIDE-3; Cell death-inducing DFFA-like effector protein C; Fat-specific protein FSP27 homolog

Immunogen:

Synthesized peptide derived from the C-terminal region of human CIDE-3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A