

Phospho-BLK (Y501) Antibody



PACO05009

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:**Gene ID:**

BLK

Uniprot

P51451

Synonyms:

BLK; Tyrosine-protein kinase Blk; B lymphocyte kinase; p55-Blk

Immunogen:

Synthesized peptide derived from human Blk around the phosphorylation site of Y501.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A