

CBLN3 Antibody



PACO05032

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:

Gene ID:

CBLN3

Uniprot

Q6UW01

Synonyms:

CBLN3; Cerebellin-3

Immunogen:

Synthesized peptide derived from the C-terminal region of human Cerebellin 3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A