

# CDH24 Antibody



PACO05033

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

ELISA:1:10000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

CDH24

**Uniprot**

Q86UP0

**Synonyms:**

CDH24; CDH11L; Cadherin-24

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human Cadherin-24.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A