

CDH22 Antibody



PACO05034

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

CDH22

Uniprot

Q9UJ99

Synonyms:

CDH22; C20orf25; Cadherin-22; Pituitary and brain cadherin; PB-cadherin

Immunogen:

Synthesized peptide derived from the Internal region of human Cadherin-22.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A