

DCLRE1C Antibody



PACO05037

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:10000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

DCLRE1C

Uniprot

Q96SD1

Synonyms:

DCLRE1C; ARTEMIS; ASCID; SCIDA; SNM1C; Protein artemis; DNA cross-link repair 1C protein; Protein A-SCID; SNM1 homolog C; hSNM1C; SNM1-like protein

Immunogen:

Synthesized peptide derived from human Artemis around the non-phosphorylation site of S516.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A