

DNAI2 Antibody



PACO05041

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

DNAI2

Uniprot

Q9GZS0

Synonyms:

DNAI2; Dynein intermediate chain 2; axonemal; Axonemal dynein intermediate chain 2

Immunogen:

Synthesized peptide derived from the N-terminal region of human Dynein IC2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A