

TNS3 Antibody



PACO05057

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

TNS3

Uniprot

Q68CZ2

Synonyms:

TNS3; TEM6; TENS1; TPP; Tensin-3; Tensin-like SH2 domain-containing protein 1; Tumor endothelial marker 6

Immunogen:

Synthesized peptide derived from the Internal region of human Tensin3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A