

RANBP17 Antibody



PACO05064

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

RANBP17

Uniprot

Q9H2T7

Synonyms:

RANBP17; Ran-binding protein 17

Immunogen:

Synthesized peptide derived from the N-terminal region of human Ran BP-17.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A