

**FOXD4/FOXD4L4/FOXD4L5/FOXD4L2/FOXD4L3/FOX
D4L6/ Antibody**
PACO05083



Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:

Gene ID:

FOXD4/FOXD4L4/FOXD4L5/FOXD4L3/FOXD4L6

Uniprot

Q12950/Q8WXT5/Q5VV16/Q6VB85/Q6VB84/Q3SYB3

Synonyms:

FOXD4; FKHL9; FOXD4A; FREAC5; Forkhead box protein D4; Forkhead-related protein FKHL9; Forkhead-related transcription factor 5; FREAC-5; Myeloid factor-alpha; FOXD4L4; FOXD4B; Forkhead box protein D4-like 4; FOXD4-like 4; Forkhead box prote

Immunogen:

Synthesized peptide derived from the Internal region of human FoxD4/D4L

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A